

Avance® Solo Adapt

NEGATIVE PRESSURE WOUND THERAPY

Avance® Solo Adapt NPWT System

Manufacturer Mölnlycke Health Care AB
Gamlestadsvägen 3C, Box 13080
SE-402 52 Göteborg, Sweden

www.molnlycke.com

www.molnlycke.com/symbols

STERILE EO	The device is sterilized using ethylene oxide
	Do not reuse
	Caution, see instructions for use
	Only use the type and model of lithium batteries specified for this product by Mölnlycke Health Care, see section 8.
	Follow instructions for use
	Temperature limitation
	Keep dry Keep away from rain
	Keep away from sunlight Keep away from heat
	Catalogue number
	Use by date / Expiry date
	Batch code
	Serial number

Rx only

Caution: Federal [US] Law restricts this device to sale or on the order of a physician (or properly licensed practitioner).

Avance Solo Adapt NPWT System is intended to be used by healthcare professionals by directions of these instructions for use. Information for patient or lay person is provided in a separate patient user manual provided by Mölnlycke Health Care. The healthcare professional shall ensure that the patient user manual is handed to the patient or lay person as appropriate.

For prescription of therapy in to home care, the prescribing healthcare professional shall confirm that the patient or lay caregiver understands how the pump and canister works and how the healthcare professional shall ensure that the patient or lay caregiver understands that the patient or lay caregiver is able to perceive audible and visual notifications and alarms from the pump and to troubleshoot trouble by guidance in the Patient User Manual. The patient or lay person shall be advised to contact prescribing healthcare professional if any concern on safe use of Avance Solo Adapt NPWT System.

1. Product description

Avance Solo Adapt Negative Pressure Wound Therapy (NPWT) System consists of Avance Solo Adapt Pump, Avance Solo Canister 50 ml, Avance Solo Adapt Foam, Avance Solo Adapt Film and Avance Solo Adapt Transfer Port which together form a system for wound management via the application of negative pressure wound therapy.

Avance Solo Adapt Pump, a battery powered single patient use pump with a 14 day lifespan, single button operated with visual and audible alarms and notifications.

Avance Solo Canister 50 ml, a single use canister attached to the pump for collection of wound fluid and exudate.

Avance Solo Adapt Foam, a single use polyurethane foam wound filler.

Avance Solo Adapt Film, a single use breathable soft silicone film, to fixate the foam and to achieve a seal.

Avance Solo Adapt Transfer Port, a single use polyolefin/polyurethane transfer port, for transporting exudate.

Avance Solo Adapt NPWT System maintains negative pressure nominally at -125 mmHg to the wound and enables exudate management by collecting exudate in Avance Solo Canister 50 ml. Avance Solo Adapt NPWT System is intended for adults.

Material content

Film: polyethylene, polyurethane, soft silicone, polycarbonate adhesive

Foam: polyurethane

Canister: polycarbonate, polyurethane

Pump: polycarbonate, acrylonitrile butadiene styrene, thermoplastic elastomer

Transfer Port and tubes with clamp: Polyurethane, polyester, polyacrylic adhesive, polyolefin based thermoplastic elastomer, polyethylene

Connectors: acrylonitrile butadiene styrene copolymer, thermoplastic olefin, polyethylene

Imprints:

- contamination

- entrapment or strangulation

- kinking or blocking of tubes

- being exposed to sources of heat

Coverage: may sharp edges or bones in the wound prior to the application of the dressing due to the risk of puncturing of organs and blood vessels.

Tendons, ligaments and nerves shall be protected from direct contact with Avance Solo Adapt Foam.

Do not place Avance Solo Adapt Foam into unexplored or blind tunnels.

If defibrillation is required, disconnect Avance Solo Adapt Pump and remove the dressing if the dressing positioning interferes.

Do not use the products on patient and/or user with known hypersensitivity to the materials/components of the products.

2. Indications for use

The Avance Solo Adapt NPWT System is indicated for patients who would benefit from wound management via the application of negative pressure wound therapy, particularly as the device may promote wound healing through the removal of exudate, infectious material.

Avance Solo Adapt NPWT System is indicated for removal of low to moderate amounts of exudate from pressure ulcers.

The healthcare professional shall ensure that the patient user manual is handed to the patient or lay person as appropriate.

For prescription of therapy in to home care, the prescribing healthcare professional shall confirm that the patient or lay caregiver understands how the pump and canister works and how the healthcare professional shall ensure that the patient or lay caregiver is able to perceive audible and visual notifications and alarms from the pump and to troubleshoot trouble by guidance in the Patient User Manual. The patient or lay person shall be advised to contact prescribing healthcare professional if any concern on safe use of Avance Solo Adapt NPWT System.

3. Contraindications

Avance Solo Adapt NPWT System is contraindicated for patients with the following conditions:

- Malignancy in the wound or margins of the wound

- Untreated and previously confirmed osteomyelitis

- Non-enteric and unexplored fistulas

- Necrotic tissue with eschar present

- Exposed nerves, arteries, veins or organs

- Exposed anastomotic site

4. Warnings

Application of NPWT may increase the risk of bleeding. If sudden or increased bleeding is observed, immediately take appropriate action to stop bleeding and seek emergency medical attention.

Patients at high risk of bleeding, such as patients receiving anticoagulant therapy or with altered hemostasis, must be monitored closely during therapy. Ensure to establish hemostasis before applying therapy.

Patients at risk of bleeding complications due to e.g. a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels, shall be monitored during therapy.

Patients undergoing NPWT need frequent supervision. Routinely check that the negative pressure therapy is active, the dressing should be re-secured and held to the wound. If it is necessary to disconnect the pump the time elapsed without negative pressure must comply with instructions from healthcare professional.

During therapy, make sure that the pump, quick connectors and tubes are positioned to eliminate the risk of:

- imprints

- contamination

- entrapment or strangulation

- kinking or blocking of tubes

- being exposed to sources of heat

Coverage: may sharp edges or bones in the wound prior to the application of the dressing due to the risk of puncturing of organs and blood vessels.

Tendons, ligaments and nerves shall be protected from direct contact with Avance Solo Adapt Foam.

Do not place Avance Solo Adapt Foam into unexplored or blind tunnels.

If defibrillation is required, disconnect Avance Solo Adapt Pump and remove the dressing if the dressing positioning interferes.

Do not use the products on patient and/or user with known hypersensitivity to the materials/components of the products.

6. Instructions for use

6.1. To consider before use

Avance Solo Adapt NPWT System is applicable for use on low to moderately exuding wounds up to the size of the Avance Solo Adapt Foam. The maximum allowed wound area shall not exceed the size of Avance Solo Adapt Foam (67 cm²/10.4 in²). Do not use more foam than what is supplied in one Avance Solo Adapt Pump Transfer Port kit.

Low exuding wounds are considered to be up to 0.6 g/cm² wound area/24 hours (0.6 g/16 in² wound area/24 hours). Moderate exuding wounds are considered to be up to 1.1 g/cm² wound area/24 hours (1.1 g/0.16 in² wound area/24 hours). 1 g exudate is considered equal to 1 ml.

13. Document the number of pieces of foam material used in the patient's notes.

Avance Solo Adapt Pump is Magnetic Resonance (MR) safe, do not take into a Magnetic Resonance Imaging (MRI) environment. Avance Solo Adapt Foam, Avance Solo Adapt Film and Avance Solo Adapt Transfer Port are MR safe. The impact of the foam, film and transfer port on Magnetic Resonance Tomography (MRT) / Magnetic Resonance Imaging (MRI) imaging artifacts is unknown.

Avance Solo Adapt Pump shall not be used in oxygen rich environments where there is a danger of explosion e.g. a hyperbaric oxygen unit or in therapy/investigations involving microwaves.

Avance Solo Adapt Pump is not suitable for use in the presence of flammable anesthetics.

For daily hygiene routines do not expose Avance Solo Adapt Pump with canister or dressing to extensive contact with water.

Avance Solo Adapt Foam and Avance Solo Adapt Film shall only be applied and changed by healthcare professional.

No modification of this device (pump, canister, tubes, foam or film) outside these instructions for use is allowed as modifications may significantly compromise the ability of the system to deliver therapy.

To minimize the risk of bradycardia, do not apply the dressing close to the vagus nerve.

Do not use oxidizing agents such as hypochlorite solutions or hydrogen peroxide prior to use of the Avance Solo Adapt Foam.

Products in Avance Solo Adapt NPWT System shall only be used with products specified in Section 1, Product description.

Do not cut the tube or detach tube from the canister.

For CT scans and X-ray investigations, keep the pump out of the X-ray or scanner range. In the event that the pump is within a CT scan or X-ray range, ensure that the pump functions correctly after the procedure.

If Avance Solo Canister or Avance Solo Adapt Pump is broken, pause the pump, replace the broken product and restart the system.

5. Precautions

Before starting therapy, evaluate the patient's nutritional status and address severe malnutrition. Also consider the patient's size and weight.

Signs of possible infection or complications must be addressed immediately. Monitor the device, wound, surrounding skin, and patient status accordingly to ensure effective and safe treatment and patient comfort.

For patients with ischemic condition or at application of a circumferential dressing, extra monitoring of wound status is required to avoid risk of compromised circulation.

Ensure that the battery lid on the Avance Solo Adapt Pump is closed during therapy.

Always clamp the transfer port tube and the canister tube before disconnecting them.

The potential for electromagnetic interference in all environments cannot be eliminated. Use caution if the pump is near electronic equipment such as anti-theft equipment or metal detectors and ensure proper function according to section 7.2 Avance Solo Adapt NPWT System - Normal use

6.4. Frequency of dressing change

Avance Solo Adapt Foam, Film and Transfer Port should be changed every 48 to 72 hours, but no less than 3 times a week, as instructed by healthcare professional.

NOTE: Ingrowth of tissue may occur if the foam is not changed around the edges to improve contact with the skin. Take care to minimize wrinkles and gaps. Consider using additional Avance Solo Adapt Film or start over and reseal, if necessary.

6.5. Dressing removal procedure

To remove the dressing, perform the following steps:

26. If the pump is active, pause therapy by pressing down the green start button, release after two (2) seconds.

27. Block both the canister tube and transfer port tube by positioning the slide clamps next to the connector and slide across the tubes until secured.

28. Disconnect the canister tube from the transfer port tube by squeezing the connector from both sides and pull apart.

29. Gently peel back one corner of the film and stretch to facilitate breakage of the seal. Proceed with this technique (in the direction of hair growth) until the film is completely removed.

30. Remove the film gently.

The film shall cover the foam and adhere to the surrounding skin with a margin of 4-5 cm/ 1.6-2.0 inches. It is preferable to apply the film as a single piece. In case several pieces are required, ensure overlap of the pieces of 4-5 cm/ 1.6-2.0 inches for proper adhesion. Minimize number of overlaps to avoid moisture build up on the film during application.

Do not place Avance Solo Adapt Foam onto intact skin or non-protected areas or retention sutures.

Do not place Avance Solo Adapt Pump with canister in water or other liquids. Disconnect the pump and canister if ingress of water is observed.

For daily hygiene routines do not expose Avance Solo Adapt Pump with canister or dressing to extensive contact with water.

Avance Solo Adapt Foam and Avance Solo Adapt Film shall only be applied and changed by healthcare professional.

No modification of this device (pump, canister, tubes, foam or film) outside these instructions for use is allowed as modifications may significantly compromise the ability of the system to deliver therapy.

To minimize the risk of bradycardia, do not apply the dressing close to the vagus nerve.

Do not use oxidizing agents such as hypochlorite solutions or hydrogen peroxide prior to use of the Avance Solo Adapt Foam.

Products in Avance Solo Adapt NPWT System shall only be used with products specified in Section 1, Product description.

Do not cut the tube or detach tube from the canister.

For CT scans and X-ray investigations, keep the pump out of the X-ray or scanner range. In the event that the pump is within a CT scan or X-ray range, ensure that the pump functions correctly after the procedure.

If Avance Solo Adapt Pump is broken, pause the pump, replace the broken product and restart the system.

6.6. Daily use

To pause therapy, press and hold the green start button, release after two (2) seconds. When paused, the pump will automatically start after 60 min.

To restart therapy, press and hold the green start button, release after two (2) seconds until the pump activates and the start button flashes.

Regularly check that negative pressure is active by monitoring visual and audible notifications and alarms from the pump. The dressing should be contracted and firm to the touch.

The tubes should be laid to prevent bends and kinks of the tubes and to avoid risk of imprints and discomfort to the patient.

7.2. Avance Solo Adapt NPWT System - Normal use

Avance Solo Adapt Pump displays the following visual and audible signals to inform the user that Avance Solo Adapt NPWT System is performing under normal operation.

AUDIBLE AND VISUAL DISPLAY	DESCRIPTION	OPERATION	COMMENT

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© Schellco Schellco USA, Inc.